Geometrisk optik

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Geometrical optics (ray optics)

Translated to English from the Swedish lab instruction. Book chapter recommended reading also translated but not usable for not Swedish capable.


Preparation reading

Read the book “Våglära och optik by Göran Jönsson”, *imaging with thin lenses* (pages 229 - 236), *the eye* (pages 271 - 273), *color and color vision* (pages 274), *additive and subtractive color mix* (pages 274 - 277), *glasses* (pages 279 - 282), *the camera* (pages 282 - 291), *angular magnification* (pages 293 - 294), *magnifying glass* (pages 294 - 295), *the microscope* (pages 295 - 297), and *binoculars* (pages 297 - 302).

Preparation questions

The solutions, or attempts to solutions must be given to the supervisor for control at the beginning of the laboratory session in order for the student to be allowed to carry out to the lab.

1. An old digital camera has a 10 mm focal length lens. The image sensor, that is the light-sensitive part of the camera, has square image elements (pixels) of size $7.5 \mu\text{m} \times 7.5 \mu\text{m}$. A test map with a variety of closely spaced line pairs (see figure to the right) is located 1.0 m from the lens. How narrow can the black bars (on the test card) be for them to be resolved by the digital camera? Note that in this case the resolution is limited only by the size of the pixels.


Answer: 0.74 mm. Convince yourself that you understand why the answer is not 0.75 mm.

2. The focal length of a positive lens can be determined as follows. If the distance L between object and image is greater than $4f$, there are two lens positions that result in a sharp image. Call the distance between the two possible positions of the lens d . (See figure on next page.) Using L and d , f can be determined. Show that

$$f = \frac{L^2 - d^2}{4 \cdot L}$$

Tip: Write down Gauss lens formula using f , a_1 and b_1 . Since $b_2 = a_1$ (for symmetry reasons, see the figure) $L = 2a_1 + d$. Also note that $b_1 = L - a_1$.

3. A positive lens is placed on an optical bench. Let the optical axis be the x -axis. The positive lens is located at $x = 0$ mm. The image created by the positive lens is then at $x = 1200$ mm. If a lens L_2 then is placed at $x = 1000$ mm, the image from L_1 will be imaged further, and the final image will be at $x = 1400$ mm. Determine the focal length of lens L_2 .

Answer: -400 mm

4. To determine the size of the image sensor and the image elements in an iPhone 5s, two crossed rulers were spaced 365 mm apart and photographed. When observing the image, 115 mm of the horizontal ruler and 86 mm of the vertical were visible. Use the table to the right to answer the questions below. The lens can be considered to be a thin lens.

- How far from the focal point of the lens should the image sensor be placed if the image is to be sharp?
- How large is the image sensor of the camera?
- The camera has square pixels. How long is the side of a pixel?

Answer: a) $47 \mu\text{m}$ b) $1.31 \text{ mm} \times 0.977 \text{ mm}$
c) $0.40 \mu\text{m}$.

Model	Focal length / mm	Number / Mpx
Apple iPhone 4	3,85	5
Apple iPhone 4S	4,28	8
Apple iPhone 5	4,10	8
Apple iPhone 5C	4,28	8
Apple iPhone 5S	4,1	8
Samsung Galaxy S3	3,7	8
Samsung Galaxy S4	4,2	13
Samsung Galaxy S5	4,8	16
Nexus 4	4,6	8
Nexus 5	4,0	8
Sony Xperia Z	4,1	13,1
Sony Xperia Z1	4,9	20,7
Sony Xperia Z2	4,9	20,7
Sony Xperia Z Ultra	3,0	8
HTC One (M8)	3,8	4
LG G2	4,0	13

Table 1. Data for some cellular phones. Thanks to Emil Ljungberg, F10.

Experiments

Tasks 1 and 2 are carried out at the long optical benches while tasks 3, 4 and 5 are carried out at the short ones.

1. Thin positive lenses

- a) Place a light source, an object, the lens labeled L1 and a screen on the optical bench. Make sure that the distance between the object and the screen is greater than 1 meter. Create a sharp image of the object on the screen by use of the equipment. Note that there are two options that can be used in order to get an image on the screen - a "projector mode" and a "camera mode". Use the projector mode and measure the distances needed to determine the focal length f_1 of L1 and the lateral magnification M of the image (the ratio of image length to object length measured perpendicular to the optical axis).
- b) Determine the focal length using the method described in Preparation question 2. Note that d is the distance between the lens projector position and its camera position.
- c) Use a method of your choice to determine the focal length and refractive strength of one of the positive eyeglass lenses in the lab.
- d) Calculate the distance $a + b$ expressed in the focal length of the lens L1 when the lateral magnification $M = -1$. Check the result using the optical bench.

2. Thin negative lenses

- a) Determine the focal length f_2 of lens L2 using lens L1 according to the method used in preparation question 3.
- b) Determine the focal length and refractive power of one of the negative eyeglass lenses in the lab.

3. The Galilean binocular

- a) Construct a Galilean binocular (theater binoculars) with an angular magnification of $G = 3$ using suitable lenses. Normalize the binoculars experimentally by calculating the binoculars length L for the lenses you have chosen, and set this distance between objective and eyepieces. Look at the wall map with the lines. The binoculars probably need to be adjusted, as the distance to the wall is short. Try to look at the lines with both eyes open (that is through the binoculars and next to it) and confirm that the angular magnification really is 3 times. Draw the rays for the image through your normalized theater binoculars in the answer booklet.
- b) Construct a binocular with an angle magnification of $G = 1/3$ using the same lenses as in the a-assignment. In what contexts is this type of binoculars used?

4. Kepler binoculars

- a) Using suitable lenses, construct a Kepler binoculars with the angular magnification $G = 4$. Normalize the binoculars by calculating the binoculars length L for the lenses you have selected, and set this distance between objective and eyepiece. Adjust the lens distances such that the lines on the wall map are sharply visible. Try to look at the line both through the binoculars and next to it and confirm that the angular magnification is 4 times. Draw the rays forming the image through your normalized Kepler binoculars in the answer booklet. What in particular distinguishes the Kepler binoculars from the Galilean binoculars?
- b) Using three positive lenses with focal lengths 50 mm, 100 mm and 200 mm, you will build a terrestrial binoculars, that is a Kepler binoculars with an image reversal lens ($M = -1$) in the

center. Start by calculating the distances to be set between lenses and make the binoculars as compact as possible!

5. The microscope

Build a normalized microscope using two positive lenses with focal lengths of 50 mm and 100 mm. The tube length (the distance between the focal points of the lens and the eyepiece) should be 160 mm. Calculate the magnification of the microscope and verify that the answer is reasonable by "using" the microscope. Use an illuminated dull glass plate as an object and adjust the microscope image by moving the object.

6. The digital camera / cellular phone

a) Determine, using the method in Preparation question 4, the size of the image sensor in a digital camera - or in your own mobile phone.

b) Calculate the size of the pixels. You can assume that they are squares.

7. Color composition of images on screen and images in print

a) In the folder "Geometrisk optik" you will find a photo called "Laos". Open the image in Paint. Use a magnifying glass to study the image on the screen. Which colors does the monitor emit? Describe the color content of a white, black and yellow area in the photo.

b) Open the "Färger" image in the program Paint. Analyze the color content of the respective fields using the pipette-tool. Use the pipette (from the tools) and click on the color you want to analyze. The three-color values R, G and B are then shown. In a fully red area, $R = 255$, $G = 0$ and $B = 0$, all red dots are lit with maximum intensity while all green and blue dots are off (minimum intensity is 0.)

What are the three colors in the middle row called, and how are they constructed?

c) Study a banknote and a printed newspaper image using a magnifying glass. Which colors are the pictures made up of? How is a banknote and newspaper printed?